

## TURN EASY SHEET (TES)



## TURN EASY SHEET (TES)

The action of turning in bed can sometimes be difficult but with the help of this new FITTED SHEET with a centre panel of low friction satin, things are different.

Getting in and out of bed is also easier.

The sides and foot end are made of good quality, easy care material to help you from slipping off the side of the bed and to give a certain amount of grip to aid movement.

For the less able, there is an optional draw sheet. It has satin on the underside and remains in place for the patient to lie on. The two low friction surfaces working together enable a CARER to turn/roll or move the patient up or down the bed. This is accomplished without handling them, thereby facilitating anti-pressure practice where regular turning is required, the long (2 metre) draw sheet can be firmly tucked-in ready for the next move. It has been found that because the system is permanently in place there is no searching for a slide sheet, or worse, doing without altogether.

Using these sheets is as compliant with tissue viability requirements as it is with moving and handling procedure. Both products are made in our factory in the UK and are fully hot washable.

situations. The ease of turning takes the strain off the carer's

hospitals and nursing homes, but also in domestic

They are found to be very useful, not only in

back.



T. 01535 682630 F:01535 682639 E:sales@thorpemill.co.uk





# Positioning aids

## Positioning aids






Used in the prone position to facilitate head, shoulder and arm movement. Made with combustion modified foam and FR wipe-clean vapour permeable cover. With carry strap. Available in 3 sizes:

## Code Length Width

	OILO	oodo	Longui	TTIGUT	Tieigitt
	Small	PWS	300 mm	300 mm	25/75 mm
	Medium	PWM	650 mm	600 mm	50/175 mm
	Large	PWL	650 mm	600 mm	50/250 mm



Postural rolls are useful in various therapeutic activities for moving weight, rolling

The T-roll positioning aid has tubular shaped CMHR foam core. With a wipe-clean

100 x 100 x 150 mm




# Code Length

		Longin	Didillotoi
Small	PRS	600 mm	200 mm
Medium	PRM	900 mm	300 mm
Large	PRL	900 mm	375 mm



T-roll positioning aid







Name





Code



Thorpe Mi

Telephone: 01535 682630 Fax 01535 682639 Email: sales@thorpemill.co.uk Web: www.thorpemill.co.uk







22 x 15 x 5 cm

Small wedge covered in vapour permeable wipe clean material.



preventing the patient from slipping down the bed at the same time reducing pressure on the heels. Covered in vapour permeable wipe clean material and has a non-slip base.

This cushion simulates the action of a knee break on a profiling bed

