




Pressure Mapping Analysis



Male 83 kg Hospital bed laid flat
Avg. Press. 19.66
Peak Press. 63.95
Min. Press. 10.34
Area (cm²) 2740.32

mmHg	
56.12	2
54.54	0
52.97	1
51.39	0
49.81	0
48.23	4
46.65	0
45.07	2
43.49	2
41.91	3
40.34	3
38.76	7
37.18	12
35.6	17
34.02	26
32.44	26
30.86	35
29.29	40
27.71	44
26.13	63
24.55	72
22.97	91
21.39	108
19.81	137
18.24	167
16.66	153
15.08	160
13.5	158
11.92	185
10.34	183
0.0	8539
10.34 - 56.12	
Sensor/Group	Value
FX100/64.160.02.S0123	
Avg Pres.	19.66
Peak Pres.	63.95
Min Pres.	10.36
Area (cm ²)	2740.32

RMS

Replacement Foam Mattress

Description

Foam mattress, castellated on both sides to aid in air circulation when in use, and to help the mattress contour to the bed frame when profiled. This mattress has a 4 way turning cycle. Two way stretch wipe clean vapour permeable P.U. cover and zip on three sides.

Dimensions

(L)1980 x (W)910 x (H)150 mm

Warranty

3 years

Max User Weight

171Kg/30 stones

Fire Retardancy

BS 7177 : 2008

1996 Medium Hazard Ignition sources 0.1.5. (Crib 5), BS 6807 : 1990