
Thorpe Mill Ltd Unit 1a Aireside Business Park, Royd Ings Ave, Keighley West Yorkshire, BD21 4BZ,

Tel: 01535 682 630 Fax: 01535 682 639

Email: sales@thorpemill.co.uk Web: www.thorpemill.co.uk





## **Optimise Community Plus profiling mattress**

4-way turn profiling mattress

Product Code: OPCOMP208514

Community Plus, our new high grade community mattress has been developed with the community hospital bed in mind. The mattress has pressure redistributing properties, which are conducive to helping to reduce the risk of pressure ulcers and has been designed to provide a high level of comfort and support at an affordable price.

The 4-way turn mattress, is comprised of a double sided castellated high-density combustion modified polyurethane, covered in a multi-stretch vapour-permeable fabric. The cover can be made with either sewn or high frequency welded seams\* and has zips on three sides, which can also be concealed with waterfall flaps if required\*.

- The high-density foam used in the Community Plus, ensures maximum rest and superior longevity.
- The castellated design is conducive to helping to reduce the risk of pressure ulcers for the patient, as it increases the pressure redistributing properties of the mattress by increasing the surface contact area, reducing shear and friction.
- 4-way turn maximises the mattresses potential and extends the life of the mattress.
- Double sided castellation, helps when profiling the bed and aids the air flow within the foam cells.
- 3-sided zip helps with washing, fitting and inspection.

\*Welded seams and waterfall flaps to conceal the zip are availiable upon request.

